


## H236 BLUE



### General Description

Blue polyurethane film on paper carrier suitable for visually detectable first aid dressings.

Property	Typical Value
Thickness (micron)	84
Weight (g/m <sup>2</sup> )	100
MVTR, upright (g/m <sup>2</sup> per 24 shrs @ 37°C/98.6°F)	2000
Tensile strength MD (gf/cm)	2000
Tensile strength TD (fg/cm)	2000
Nominal length at max width (lm)	1000
Maximum width (mm)	1300
Minimum width (mm)	50
Carrier	92 gsm double sided release paper
Liner	N/A
Adhesive	N/A
Stripping load (release) (gf/cm)	1.0
Elongation at break, MD (%)	900
Elongation at break, TD (%)	900