

PBA74/A16 BEIGE



General Description

Beige medical grade polyurethane/urethane blend, breathable embossed film with A16 adhesive.

Property	Typical Value
Thickness (micron)	152
Film weight (g/m ²)	74
MVTR, upright (g/m ² per 24 hrs @ 37°C/98.6°F)	500
Nominal length at max width (lm)	500
Maximum width (mm)	1190
Minimum width (mm)	40
Carrier	N/A
Liner	81 gsm single sided release glassine paper
Adhesive	A16
Stripping load (release) (gf/cm)	5.0
Other information	Minimum adhesion to steel 150 gf/cm