

## PBA74 BLUE



### General Description

Blue, medical grade polyurethane/urethane blend, suitable for visually detectable first aid dressings.

Property	Typical Value
Thickness (micron)	152
Film weight (g/m <sup>2</sup> )	74
MVTR, upright (g/m <sup>2</sup> per 24 shrs @ 37°C/98.6°F)	1600
Tensile strength MD (gf/cm)	1000
Tensile strength TD (gf/cm)	500
Nominal length at max width (1m)	1500
Maximum width (mm)	1220
Minimum width (mm)	250
Carrier	Net interleave only
Liner	N/A
Adhesive	N/A
Elongation at break, MD (%)	240
Elongation at break, TD (%)	300